


WaffleGrip Trendelenburg Positioning System (PN A300)

a HotDog® Patient Warming accessory

Instructions for Use



WaffleGrip
Warm & Secure



A300 (10 Kits)

Each Kit Contains:

- (1) Waffle Grip Pad
- (1) Draw Sheet (Patient Lifting and Arm-Tucking)
- (1) Safety strap

Kit does not include:

- U300 Warming Mattress
- (6) Side rail securing straps

U300 HotDog Warming Mattress

Operating Table Mattress

Operating Table

Cross-Section

Reinforced Plastic Flap --Tuck Between HotDog Warming Mattress and Operating Table Mattress

INSTRUCTIONS FOR USE

1. Place U300 HotDog Patient Warming Mattress on operating table with sensor and cable toward head end (a). Align the perineal cutout on the warming mattress with that of the operating table (b). Attach mattress to side rails using **at least 4** straps as shown (c).
2. Place non-slip WaffleGrip over the HotDog warming mattress as shown in label, and unfold. Feed mattress cable through pad porthole (a). Align sealed fold lines on WaffleGrip with ends of warming mattress, with reinforced plastic flap at perineal cutout.
IMPORTANT -- TUCK REINFORCED PLASTIC FLAP BETWEEN WARMING MATTRESS AND OPERATING TABLE MATTRESS [SEE CROSS-SECTION (b)]. Flap should not be visible or hanging amongst perineal cutout. Tuck remaining WaffleGrip pad (beyond fold lines) under warming mattress.
3. Place Non-Woven Draw Sheet transverse across table as shown, with textured WaffleGrip toward patient (a).
4. Connect HotDog patient warming mattress to WC52 or WC77 HotDog Controller using A112 blue cable. Select temperature. See Warming Mattress IFU and Controller IFU for instructions on operating HotDog patient warming mattress.
5. **IMPORTANT** -- PATIENT MUST BE IN DIRECT SKIN CONTACT WITH WAFFLEGRIP to prevent slipping in Trendelenburg. Position patient using draw sheet. Align patient's buttocks over tucked end of WaffleGrip. Tuck arms as shown or per hospital positioning protocol (a).
6. Secure Safety Strap across patient's upper chest by looping straps under side rails and feeding through metal brackets on each side. Secure straps over patient's chest (a). Strap can be used to help secure HotDog warming blanket.
7. See separate instructions for most effective HotDog Patient Warming Blanket positioning during Trendelenburg cases.

USE ONLY AS INDICATED

ATTENTION: Waffle pattern will be temporarily visible on skin, but will dissipate shortly following procedure.

IMPORTANT:

1. Do not use lubricants, oils or powders on patient's skin in contact with WaffleGrip. Wipe skin with alcohol if any exist. (A decrease in friction will reduce the non-slip capability of the pad)
2. Not designed for use in Reverse Trendelenburg.



Augustine Surgical Inc.
Making Surgery Safer



Augustine Temperature Management, LLC
Eden Prairie, MN 55344 USA
www.hotdogwarming.com Ph: 952-465-3500

