



air-free
PATIENT WARMING

Introducing the Most Advanced & Versatile Patient Warming Blanket

Air-Free | Better Warming | Cost Effective B500 Universal Blanket

NEW & IMPROVED



Versatile. Efficient. Effective.



Warm the legs, torso, arms, or whatever is necessary to achieve normothermia.



- A typical operating room will need two (2) B500 Universal Warming Blankets in order to successfully warm all procedures.
- A WC77 HotDog Temperature Management Controller can operate one mattress and four blankets simultaneously.
- Combine with a warming Mattress for best outcomes.
- Can be connected or overlapped for versatile positioning.



Augustine Surgical Inc.
Making Surgery Safer

www.hotdogwarming.com
©2022 Augustine Temperature Management. All Rights Reserved.
Patents: hotdogwarming.com/patents
Rx Only
M128EN Rev B 6/22

6581 City West Parkway
Eden Prairie, MN 55344
952-746-1720 or 888-439-2767
cs@aug surg.com